



Low Volume dialysis cell

6.2729.200

Including four 6.2744.060 stoppers. For rapid Inline Dialysis.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.



Scope of delivery 6.2729.200

Qt.	Order no.	Description
4 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optional accessories

Order no.	Description	
6.2714.010	Dialysis membrane (cellulose acetate) Nominal pore size : 0.2 µm. Set of 50 pieces. For inline dialysis	
6.2714.030	Dialysis membrane (polyamide) For inline dialysis and filtration	
6.2744.000	Pressure screw PVDF With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)	