



## EMPTY sensor for level sensor for 1 L bottles

6.2769.100

This sensor, mounted in the 6.1626.000 bottle attachment, results in the level sensor EMPTY for 1 L bottles. Included are two sensor rods, one aspiration tubing and two holders for fixing the tubing to one of the sensor rods.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Scope of delivery 6.2769.100

Qt.	Order no.	Description
1 PCS	6.1829.010	FEP aspiration tube M6, 0.25 m With M 6 thread. For Dosing Unit

