



Set for exact temperature adjustment

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Scope of delivery 6.5616.100

Metrohm AG

1 PCS

6.1111.010

Pt100 temperature sensor for stability measuring devices

Pt100 temperature sensor (length 17.5 cm) with fixed cable and mini-DIN-plug. This glass-free temperature sensor is well suited for use with a 711 Liquino and instruments for stability measurement (743 Rancimat, 763 PVC Thermomat, 892 Professional Rancimat, or 895 Professional PVC Thermomat).



2 PCS

6.2042.040

Clamp for temperature sensor

Holding clamp for Pt100 temperature sensor for stability measuring instruments



1 PCS

6.2326.000

Silicone oil for stability measuring instruments (50 mL)

Silicone oil for the determination of temperature correction in stability measuring instruments



1 PCS

6.2418.000

Air tube for oil measurements

Set of 12 pieces. In combination with 6.1429.040.

