



# "ProfIC Vario 6 Anion" – Professional IC Vario system with Inline Dilution and Inline Ultrafiltration

## ProfIC Vario 6 Anion

The Professional IC Vario system with Metrohm Inline Dilution, Metrohm Inline Ultrafiltration, sequential suppression and conductivity detection enables fully automatic determination of anions or cations that are present in very high concentrations and that have matrices containing particles.


Typical areas of application:

- Rinse, process, and waste water
- Extractions and digestion solutions
- Food samples

Schematic representation

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Scope of delivery ProfIC Vario 6 Anion

| Qt.   | Order no.  | Description                                                                                                                                                                                                                                                                                                                                                                                                                    |                                                                                       |
|-------|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 1 PCS | 2.741.0010 | 741 Magnetic Stirrer                                                                                                                                                                                                                                                                                                                                                                                                           |                                                                                       |
|       |            | Magnetic stirrer for sample changers.                                                                                                                                                                                                                                                                                                                                                                                          |                                                                                       |
|       |            | The magnetic stirrer is located beneath the sample rack.                                                                                                                                                                                                                                                                                                                                                                       |    |
| 1 PCS | 2.800.0010 | 800 Dosino                                                                                                                                                                                                                                                                                                                                                                                                                     |                                                                                       |
|       |            | The 800 Dosino is a drive with write/read hardware for intelligent Dosing Units. With fixed cable (length 150 cm).                                                                                                                                                                                                                                                                                                             |   |
| 1 PCS | 2.850.9010 | IC Conductivity Detector                                                                                                                                                                                                                                                                                                                                                                                                       |                                                                                       |
|       |            | Compact and intelligent high performance conductivity detector for intelligent IC instruments. Outstanding temperature stability, the complete signal processing within the protected detector block and the latest generation of DSP – Digital Signal Processing – guarantee the highest precision of the measurement. No change of measuring ranges (not even automatic ones) is required, due to the dynamic working range. |  |
| 1 PCS | 2.858.0020 | 858 Professional Sample Processor – Pump                                                                                                                                                                                                                                                                                                                                                                                       |                                                                                       |
|       |            | The 858 Professional Sample Processor – Pump processes samples from 500 µL to 500 mL. The sample transfer takes place either with the installed bidirectional two-channel peristaltic pump or with an 800 Dosino.                                                                                                                                                                                                              |  |

---

1 PCS

2.940.1500

940 Professional IC Vario ONE/SeS/PP

The 940 Professional IC Vario ONE/SeS/PP is the intelligent IC instrument with sequential suppression and a peristaltic pump for suppressor regeneration. The instrument can be used with any separation and detection methods.

Typical areas of application:

- Anion or cation determinations with sequential suppression and conductivity detection
- Trace analysis for anions or cations
- Online monitoring for anions or cations



---


1 PCS

6.05330.110

IC equipment: Inline Ultrafiltration 2 - push mode

Accessory set for Inline Ultrafiltration 2 - push mode. For use with the 858 Professional Sample Processor/ 919 IC Autosampler plus.

Recommended set for the combination with Inline Dilution (MIDT).



---

1 PCS

6.5330.120

IC equipment: Inline dilution

Accessory set for Inline Dilution. For use with the 858 Professional Sample Processor, 800 Dosino and 741 Magnetic Stirrer.



## Optional accessories

| Order no.  | Description                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|------------|--------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 2.941.0010 | 941 Eluent Production Module         | <p>The 941 Eluent Production Module enables the automatic production of an eluent. It enables continuous working without manual intervention and guarantees stable retention times. It can be combined with all Metrohm IC instruments and is monitored by MagIC Net.</p>                                                                                                                                                                             |
| 6.2041.440 | Sample rack 148 x 11 mL + 3 x 300 mL | <p>Sample rack for 148 x 11 mL samples with 3 rinsing beakers.</p>                                                                                                                                                                                                                                                                                                                                                                                    |
| 6.2041.760 | Sample rack 54 x 11 mL + 1 x 300 mL  | <p>Sample rack with small diameter for 54 x 11 mL sample vessels in two rows and additionally 1 x 300 mL bottle.</p>                                                                                                                                                                                                                                                                                                                                |
| 6.2743.050 | Sample vessel 11 mL                  | <p>2000 pieces. For IC Sample Processors and VA Autosamplers.</p> <p>The IC sample vessels are made of polypropylene (PP) plastic for health-care applications. Production and packaging both take place under sterile and clean room conditions (Clean Room Class 7, in accordance with DIN ISO 14644). This guarantees minimal degrees of contamination. The sample vessels are subjected to random checks for leachable anions and cations.</p>  |
| 6.2743.070 | Stopper with perforation             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |



---

6.2744.080 M6 thread / UNF 10/32 adapter

With M6 outer thread. For IC instruments.



---

6.2832.000 MSM Rotor A

Suppressor rotor for all IC instruments with MSM (Metrohm Suppressor Module)



---

6.2842.000 MSM-HC Rotor A

Suppressor rotor for all IC instruments with MSM-HC (Metrohm Suppressor Module with high capacity)



---

6.2842.020 Adapter sleeve, MSM

Not only the MSM Rotor A but also the MSM-LC Rotor A can be installed in the rotor housing with the aid of the adapter sleeve.



---

6.2842.200

MSM-HC Rotor C

Suppressor rotor for all IC instruments with MSM-HC drive (Metrohm Suppressor Module with high capacity) for cation suppression.



---

6.2844.000

MSM-LC Rotor A

Suppressor rotor for use with 2 mm columns.



---

6.6059.422

MagIC Net 4.2 Professional: 1 license

Professional chromatography software from Metrohm for controlling all intelligent Professional IC systems, Compact IC systems and their peripherals, all detectors and the various Auto samplers, 800 Dosino, 771 Compact Interface, etc. The software permits checks, data acquisition, data evaluation and data monitoring as well as report generation of ion chromatographic analyses.



Graphical user interface for routine operations, extensive database programs, method development, configuration and manual control panel; very flexible user management, efficient database operations, extensive data export functions, individually configurable report generator, control and monitoring of all system components and the chromatography results.

MagIC Net Professional complies fully with FDA Regulation 21 CFR Part 11 as well as GLP.

MagIC Net is available in 15 dialog languages: German, English, Chinese, Traditional Chinese, French, Italian, Spanish, Portuguese, Bulgarian, Czech, Hungarian, Japanese, Korean, Russian, Slovakian

With this version, the graphic user interface has been modernized, which improves, among other things, the scaling and sharpness of the user interface on 4K monitors.

1 license

6.6059.423

#### MagIC Net 4.2 Multi: 3 licenses

Client/server-compatible chromatography software from Metrohm for controlling all intelligent Professional IC systems, Compact IC systems and their peripherals, all detectors and the various Autosamplers, 800 Dosino, 771 Compact Interface, etc.

The software permits checks, data acquisition, data evaluation and data monitoring as well as report generation of ion chromatographic analyses. Data and method are centrally stored on a Server.



Graphical user interface for routine operations, extensive database programs, method development, configuration and manual control panel; very flexible user management, efficient database operations, extensive data export functions, individually configurable report generator, control and monitoring of all system components and the chromatography results.

MagIC Net Multi complies fully with FDA Regulation 21 CFR Part 11 as well as with GLP.

MagIC Net is available in 15 dialog languages: German, English, Chinese, Traditional Chinese, French, Italian, Spanish, Portuguese, Bulgarian, Czech, Hungarian, Japanese, Korean, Russian, Slovakian

With this version, the graphic user interface has been modernized, which improves, among other things, the scaling and sharpness of the user interface on 4K monitors.

Performance for client/server systems with servers outside the in-house company network has been improved.

Client/server version with 3 licenses

---