



μStat-i 400 Bipotentiostat/Galvanostat /Impedance Analyzer (EIS)

DRP-STAT-I-400

μStat-i 400 is a portable BiPotentiostat/Galvanostat/Impedance Analyzer (EIS) that can be applied for Voltammetric, Amperometric, Potentiometric and EIS measurements and can be used with one- or two-working electrodes configuration. Li-ion Battery powered, it can be easily connected to a PC via USB and wireless.

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery DRP-STAT-I-400

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

DRP-I-CABSTAT

μStat Cable connector for μStat-i 400 and
conventional electrodes

μStat Cable connector for μStat-i 400 and conventional electrodes

