



Módulo de bombeo peristáltico (2 canales)

2.1016.0010

Módulo para el montaje en el soporte de módulo del Pick&Place del OMNIS Sample Robot. Esta puesto de trabajo está compuesta por una bomba de lavado y una de succión. Se usan para purgar los sensores de un módulo Pick&Place y para vaciar el vaso de muestra después del análisis, antes de que vuelva a ser colocado en la gradilla.

Partes/accesorios 2.1016.0010

Qt.	Order no.	Descripción
1 PCS	6.01543.100	Aspiration tip automation

Aspiration tip for automated OMNIS titrations.



1 PCS

6.01803.000

Silicone tubing 8/12 2.5 m

Silicone tubing for connecting to the drainage of the Pick&Place module or pump module in order to channel any solution overflows that may occur safely into the waste canister.



1 PCS

6.01804.500

Tubing adapter for drainage tubing

Tubing adapter for drainage tubing on Pick&Place modules and pump modules in the OMNIS Sample Robot Pick&Place.




3 PCS

6.01805.000

FEP tubing connection / M6 / 66.5 cm

With light protection and protection against kinks.



1 PCS

6.01805.010

PTFE tubing connection / M8 / 66.5 cm

With light protection and protection against kinks.



1 PCS

6.05001.010

"Peristaltic" (2-channel) pump module accessory set

The accessory set contains all of the components required for the tube connections between a rinsing pump and an aspiration pump at the "Peristaltic" pump module. In addition to preassembled tubing with Luer connectors, each of these also includes two press clamps and pump tubing.



2 PCS

6.1621.000

Container 10 L

As rinsing or waste container in automated systems.



2 PCS

6.1621.100

Cap for 10-L container

Cap for 10-L container. For easy and safe transportation of filled containers.



2 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



1 PCS

6.2709.070

Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders



3 PCS

6.2740.020

Spray nozzle

For the fine-spraying of the rinsing solution.



Accesorios opcionales

Order no.	Descripción	
6.05700.010	Consumable Kit OMNIS Pump Module (2-channel)	
	Two pump tubings each and associated press clamps for replacement in the OMNIS pump module. One Consumable Kit can be used for re-equipping two peristaltic pumps of the pump module.	
6.05800.030	Upgrade peristaltic pump module	
	This set contains all of the components that are required for the upgrade of a peristaltic (2-channel) pump module to a peristaltic (4-channel) pump module in the OMNIS Sample Robot Pick&Place. In addition to the pumps required, all the necessary tubing and canisters are included to equip one additional rinsing pump and one additional aspiration pump each.	
6.05800.080	Upgrade peristaltic pump module	
	This set contains all of the components that are required for the upgrade of a peristaltic (2-channel) pump module to a peristaltic (4-channel) pump module in the OMNIS Sample Robot Pick&Place. In addition to the pumps required, all the necessary tubing and canisters are included to equip one additional rinsing pump and one additional aspiration pump each. Compatible with OMNIS software version 2.13.2 and 2.17 or higher.	