



823 Membrane Pump Unit

2.823.0010

Bomba de membrana para conexión directa a los procesadores de muestras.

Partes/accesorios 2.823.0010

Qt.	Order no.	Descripción
1 PCS	1.823.0010	823 Membrane Pump Unit
		Membrane pump for direct connection to Sample Processors.



2 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



1 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.

