


824 Easy Sample Changer

2.824.0020

El cambiador de muestras económico para aplicaciones de rutina. Easy Sample Changer Micro: cambiador de muestras para valoraciones con pequeños volúmenes de muestra. Se pueden procesar en una pasada hasta 24 muestras con un volumen máximo de 75 mL. Con el equipo KFT, el Easy Sample Changer se convierte en un cambiador de muestras Karl Fischer.

Partes/accesorios 2.824.0020

Qt.	Order no.	Descripción
1 PCS	1.741.0010	741 Magnetic Stirrer
		Magnetic stirrer for sample changers.
		The magnetic stirrer is located beneath the sample rack.
		
24 PCS	6.1432.210	Sample beaker, clear glass, 75 mL, 1 piece
		Sample beaker, clear glass, 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, 6.2041.800 and 6.2041.900, 1 piece
		
1 PCS	6.1453.250	Sample beaker, plastic (PP), 250 mL, 1 piece
		Sample beaker, plastic (polypropylene, PP), 250 mL, suitable for sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020, and 6.2068.030, 1 piece
		

1 PCS

6.1458.020

Titration head / 4 x M10 connectors

With 4 M10 screw nipples, including stoppers.



1 PCS

6.1805.110

FEP tubing / M6 / 80 cm

With light and kink protection



24 PCS

6.1903.030

Stirring bar / 25 mm

Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 PCS

6.2034.020

Mounting bracket for Magnetic Stirrer 741

To attach the magnetic stirrer to the «824 Easy Sample Changer



1 PCS

6.2041.340

Sample rack, 24 x 75 mL, (PVC)

Sample rack for 810 Sample Processor, 814 USB Sample Processor, 855 Robotic Titrosampler "Light", suitable for 24 sample beakers. The following sample beakers can be used: 6.1432.210, 6.1459.400.

Plastic: Polyvinyl chloride (PVC)



1 PCS

6.2621.140

Hexagon key 2.5 mm




1 PCS

6.2751.080

Safety shield

Safety shield for single-tower Sample Processors without 786 Swing Head.



1 PCS

6.2752.010

Plug cover

Used with the 717, 730, 760 Sample Changers, IC Sample Processors, 774 Oven Sample Processor



Accesorios opcionales

Order no.	Descripción
6.1453.220	Sample beaker, plastic (PP), 200 mL, 1 piece Sample beaker, plastic (polypropylene, PP), 200 mL, suitable for sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020, and 6.2068.030, 1 piece

