

pH Lab Irda 240V AUS

2.827.0017

Partes/accesorios 2.827.0017

Qt.	Order no.	Descripción
1 PCS	6.2013.010	Clamping ring
For support rods with a diameter of 10 mm.		
		
1 PCS	6.2016.050	Support rod / 300 mm
		
1 PCS	6.2021.020	Electrode holder
Electrode holder for 4 electrodes and 2 buret tips		
		

1 PCS

6.2161.040

Plug-in power supply unit AUS 240V/6V DC

Power supply for 827 pH lab



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors

