



## Electrodo Gel con Pt1000

6.00221.600

Electrodo pH (electrolito gel) sin mantenimiento con sensor de temperatura (Pt1000) integrado para medidas de pH rutinarias en muestras similares. Este electrodo se guarda en una solución con  $c(\text{KCl}) = \text{sat.}$  (6.2308.000) y no es apto para soluciones pobres en iones.

El indicador de envejecimiento muestra la necesidad de reemplazar el electrodo en una etapa temprana.

### Partes/accesorios 6.00221.600



Qt.	Order no.	Descripción
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm, polyethylene.



Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



## Accesorios opcionales

Order no.	Descripción	
6.2104.600	Electrode cable for plug in head U/plug F, 2x2 mm B, 1m  For connecting electrodes with Metrohm plug-in head U to Metrohm instruments (socket F).	
6.2104.610	Electrode cable for plug in head U/plug F, 2x2 mm B, 2m  For connecting electrodes with Metrohm plug-in head U to Metrohm instruments (socket F).	
6.2308.000	Electrolyte KCl (sat.), 250 mL  Electrolyte KCl (sat.), 250 mL	