



Célula de medida de la conductividad de 5 anillos con $c = 0,7 \text{ cm}^{-1}$ con Pt1000 (cable fijo de 0,65 m)

6.00925.100

Célula de medida de la conductividad de 5 anillos con constante de célula $c = 0,7 \text{ cm}^{-1}$ (valor guía), con sensor de temperatura Pt1000 integrado y cable fijo (0,65 m) para su conexión al OMNIS Measuring Module Conductivity.

Este sensor es apto para medidas de conductividades intermedias (desde $5 \mu\text{S}/\text{cm}$ hasta $20 \text{ mS}/\text{cm}$), por ejemplo, en las siguientes sustancias:

- Agua potable
- Aguas superficiales
- Aguas residuales

Partes/accesorios 6.00925.100

Qt.	Order no.	Descripción
-----	-----------	-------------

1 PCS


6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



Accesorios opcionales

Order no.	Descripción
6.02101.020	Measuring module conductivity Measurement channel for one OMNIS Titrator or Titration Module for the connection of conductivity measuring cells. 
6.2301.060	KCl ion standard 250 mL KCl ion standard ($c(\text{ion}) = 0.1000 \pm 0.0005 \text{ mol/L}$), can also be used as standard solution for conductometry 12.87 mS/cm (25°C). 
6.2324.010	Conductivity standard $100 \mu\text{S/cm}$ 250 mL Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$. 
6.2324.110	Conductivity standard $100 \mu\text{S/cm}$, 5 x 30 mL Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$. 