



Metrosep A Supp 17 - 100/4,0

6.01032.410

Las columnas Metrosep A Supp 17 son columnas de separación de aniones para su uso a temperatura ambiente con muy buena relación calidad-precio. La Metrosep A Supp 17 - 100/4,0 permite la separación de los aniones estándar en poco tiempo. Gracias a la alta flexibilidad en el flujo (hasta 1,8 mL/min) se pueden conseguir, con dependencia del problema de separación, tiempos de análisis muy cortos.

Partes/accesorios 6.01032.410

Qt.	Order no.	Descripción
2 PCS	6.2744.060	Threaded stopper
For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.		



Accesorios opcionales

Order no.	Descripción	
6.01032.500	Metrosep A Supp 17 Guard/4.0	
	The Metrosep A Supp 17 Guard/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily onto the analytical column. No tools are required.	
6.01032.510	Metrosep A Supp 17 S-Guard/4.0	
	The Metrosep A Supp 17 S-Guard/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. The guard column is simply connected to the separation column using capillary connections. No tools are required.	
6.01032.530	Metrosep A Supp 17 S-Guard - 50/4.0	
	The Metrosep A Supp 17 S-Guard 50/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. The guard column is simply connected to the separation column using capillary connections. No tools are required.	