



Metrosep A Supp 17 - 150/4,0

6.01032.420

La columna de separación Metrosep A Supp 17 - 150/4,0 es la columna predilecta para las determinaciones de aniones que precisan una buena eficacia de separación y tiempos de separación cortos a temperatura ambiente. El máximo flujo de 1,4 mL/min proporciona en estos casos la posibilidad de optimizar la determinación. Las columnas Metrosep A Supp 17 convencer por su buena relación calidad-precio.

Partes/accesorios 6.01032.420

Qt.	Order no.	Descripción
2 PCS	6.2744.060	Threaded stopper
For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.		



Accesorios opcionales

Order no.	Descripción	
6.01032.500	Metrosep A Supp 17 Guard/4.0	
	The Metrosep A Supp 17 Guard/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily onto the analytical column. No tools are required.	
6.01032.510	Metrosep A Supp 17 S-Guard/4.0	
	The Metrosep A Supp 17 S-Guard/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. The guard column is simply connected to the separation column using capillary connections. No tools are required.	
6.01032.530	Metrosep A Supp 17 S-Guard - 50/4.0	
	The Metrosep A Supp 17 S-Guard 50/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. The guard column is simply connected to the separation column using capillary connections. No tools are required.	