



Metrosep A Supp 18 - 250/4,0

6.01033.430

La selectividad de las columnas Metrosep A Supp 18 está específicamente diseñada para trabajar con eluyentes de hidróxido. La Metrosep A Supp 18 - 250/4,0 es la versión más larga de columna A Supp 18. Gracias al pequeño tamaño de las partículas, de 3,5 μm , se consiguen altas eficiencias de separación. Los picos simétricos permiten determinar fácilmente los aniones, incluso en gamas de concentración bajas de $\mu\text{g/L}$. La Metrosep A Supp 18 - 250/4,0 es especialmente adecuada para la separación de los oxohalogenuros clorito, bromato y clorato además de los aniones estándar a temperatura ambiente.

Accesorios opcionales

Order no.	Descripción
-----------	-------------

The Metrosep A Supp 18 Guard/4.0 reliably protects the Metrosep A Supp 18 analytical separation column against contamination from sample or eluent. It contains the same separation material as the Metrosep A Supp 18 IC Columns, is also made of PEEK, and is screwed directly onto the respective separating column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, without influencing chromatographic separating efficiency.

