



Metrosep C 6 - 150/2,0

6.01051.220

El material C 6 tiene una gran capacidad y convierte a la columna Metrosep C 6 - 150/4,0 versión microbore en la solución óptima para separar cationes estándar con grandes diferencias de concentración en tiempos de retención razonables. Las aguas potables con bajos contenidos de amonio pueden determinarse con esta columna.

La columna es apta para el uso en el acoplamiento IC-MS.

Partes/accesorios 6.01051.220

Qt.	Order no.	Descripción
2 PCS	6.2744.060	Threaded stopper
For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.		



Accesorios opcionales

Order no.	Descripción
6.01051.600	Metrosep C 6 Guard/2.0 The Metrosep C 6 Guard/2.0 contains the C 6 column material and is used to protect against particles and contamination. This considerably lengthens the service life of the analytical separation column. The Metrosep C 6 Guard/2.0 works according to the "On Column Guard System" and is screwed directly onto the respective separation column with virtually no dead volume.

