



# Metrosep C Supp 2 - 100/4,0

6.01053.410

El material de separación de la Metrosep C Supp 2 está basado en un copolímero de divinilbenceno-poliestireno con grupos carboxilo. Es adecuada para la determinación y separación de cationes monovalentes y bivalentes. La Metrosep C Supp 2 - 100/4,0 es la columna de separación más corta de la familia Metrosep C Supp 2. Es particularmente adecuada para la determinación de cationes estándar. Es posible alcanzar límites de determinación muy bajos en la gama de  $\mu\text{g} / \text{L}$  gracias al ruido de línea de base posterior a la supresión secuencial extremadamente bajo.

## Partes/accesorios 6.01053.410

Qt.	Order no.	Descripción
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



## Accesorios opcionales

Order no.	Descripción
6.01053.500	Metrosep C Supp 2 Guard/4.0

The Metrosep C Supp 2 Guard/4.0 contains the C Supp 2 column material and is used to protect Metrosep C Supp 2 cation columns. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The Metrosep C Supp 2 Guard/4.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.

