




Gradilla de muestras OMNIS, 16 x 120 mL, (PP)

6.02041.030

Gradilla de muestras OMNIS para la función Pick&Place del OMNIS Sample Robot apta para 16 vasos de muestra. Se pueden utilizar los siguientes vasos de muestra: 6.01400.200, 6.01400.300, 6.01400.303.

Plástico: polipropileno (PP)

Partes/accesorios 6.02041.030

Qt.	Order no.	Descripción
1 PCS	6.02042.000	Holder for rack label (3x)

The holder for rack labels permits durable labeling of the rack with an unambiguous name.



Accesorios opcionales

Order no.	Descripción	
6.01400.200	Sample beakers, plastic (PP), 120 mL, 20 pieces Sample beakers, plastic (polypropylene, PP), 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.	
6.01400.300	Sample beakers, clear glass, 120 mL, 20 pieces Sample beakers, clear glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces.	
6.01400.303	Sample beakers, amber glass, 120 mL, 20 pieces Sample beakers, amber glass, 120 mL, suitable for 602041030 OMNIS sample rack, 20 pieces	
6.01404.030	Beaker adapter for OMNIS 120 mL sample beaker Beaker adapter for centering the 120 mL sample beaker in the slide of the Pick&Place modules	