



Electrodo de doble hilo de platino para coulometría

6.0341.100

Electrodo indicador que se utiliza para la titulación Karl Fischer coulométrica.

Partes/accesorios 6.0341.100

Qt.	Order no.	Descripción
1 PCS	6.2008.040	Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Accesorios opcionales

Order no.	Descripción	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	