



Electrodo de membrana de polímero combinado, Ca

6.0510.100

Electrodo selectivo de calcio combinado con membrana de polímero.

Este electrodo ion-selectivo es adecuado para:

- Medidas de iones de Ca^{2+} ($5 \cdot 10^{-7}$ hasta 1 mol/L) en soluciones acuosas
- Titulaciones (inversas) complexométricas (por ejemplo, determinación de la dureza del agua)

Gracias a un eje de plástico robusto/irrompible de polipropileno y la protección contra impactos de la membrana de polímero, este sensor puede soportar altas cargas mecánicas.

Como electrolito de referencia se utiliza $c(\text{NH}_4\text{NO}_3) = 1 \text{ mol/L}$.

Partes/accesorios 6.0510.100

Qt.	Order no.	Descripción
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Accesorios opcionales

Order no.	Descripción	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
6.2327.000	Electrolyte NH4NO3 1 mol/L (50 mL) Electrolyte solution c(NH4NO3) = 1 mol/L (for combined Ca-ISE 6.0510.100)	