




Kit de consumibles de juntas para tapa KF Dis-Cover 6.02710.010

6.05700.170

Juntas para repuesto en tapa 6.02710.010, 16 unidades

Accesorios opcionales

Order no.	Descripción
6.02710.010	KF Dis-Cover lid for OMNIS 120 mL sample beaker, 16 pieces Lid with seals for volumetric Karl Fischer Titration in 6.01400.200, 6.01400.300 and 6.01400.303 sample beakers in the OMNIS Sample Robot Pick&Place system.

