



Electrodo de referencia LL-Ag/AgCl para CVS

6.0728.130

Sistema de referencia interno con electrolito $c(\text{KCl}) = 3 \text{ mol/L}$. Para el uso con un recipiente de electrolito de plástico 6.1245.010.

Accesorios opcionales

Order no.	Descripción
6.1245.010	Electrolyte vessel made of plastic
	Electrolyte vessel made of plastic with ceramic diaphragm for double junction reference electrode.



6.2112.130

Adapter plug contact 5.5 mm / screw contact, gold-plated

For VA electrodes

