




Célula de medida de la conductividad con $c = 0,1 \text{ cm}^{-1}$ con Pt1000 (cable fijo)

6.0914.040

Célula de medida de la conductividad de acero fino con constante de célula $c = 0,1 \text{ cm}^{-1}$, con sensor de temperatura Pt1000 integrado y con cable fijo (4 x B (4 mm)).

Este sensor es adecuado para medidas de conductividades muy bajas (desde 0 hasta $300 \mu\text{S}/\text{cm}$), por ejemplo, en agua desionizada o para medidas según la USP <645>.

El sensor se puede utilizar con la caja de adaptadores (6.2103.160) del 856 Conductivity Module.

Partes/accesorios 6.0914.040

Qt.	Order no.	Descripción
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



Accesorios opcionales

Order no.	Descripción
-----------	-------------

6.1420.100 Flow-through measuring vessel
For 12 mm dia. electrodes (pH, metal, conductivity)



6.2324.010 Conductivity standard 100 μ S/cm 250 mL
Conductivity standard for calibration of conductivity measuring cells with cell constant = 0.1/cm.



6.2324.110 Conductivity standard 100 μ S/cm, 5 x 30 mL
Conductivity standard for calibration of conductivity measuring cells with cell constant = 0.1/cm.

