



Célula de medida de la conductividad de 5 anillos con $c = 0,7 \text{ cm}^{-1}$ con Pt1000 (cable fijo)

6.0915.100

Célula de medida de la conductividad de 5 anillos con constante de célula $c = 0,7 \text{ cm}^{-1}$ (valor guía), con sensor de temperatura Pt1000 integrado y cable fijo (1,2 m) para su conexión al 856 Conductivity Module.

Este sensor es apto para medidas de conductividades intermedias (desde $5 \mu\text{S}/\text{cm}$ hasta $20 \text{ mS}/\text{cm}$), por ejemplo, en:

- Agua potable
- Aguas superficiales
- Aguas residuales

Partes/accesorios 6.0915.100

Qt.	Order no.	Descripción
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Accesorios opcionales

Order no.	Descripción	
6.2301.060	KCl ion standard 250 mL	
	KCl ion standard ($c(\text{ion}) = 0.1000 \pm 0.0005 \text{ mol/L}$), can also be used as standard solution for conductometry 12.87 mS/cm (25°C).	
6.2324.010	Conductivity standard $100 \mu\text{S/cm}$ 250 mL	
	Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$.	
6.2324.110	Conductivity standard $100 \mu\text{S/cm}$, 5 x 30 mL	
	Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$.	
6.2763.000	Flow-through cell for the 6.0915.100/6.0920.100	
	Flow-through cell for the 6.0915.100/6.0920.100 5-ring conductivity measuring cell	