



Celda de medida de la conductividad con $c = 0,5 \text{ cm}^{-1}$ con Pt1000 (cable fijo)

6.0917.080

Celda de medida de la conductividad de 4 hilos con constante de celda $c = 0,5 \text{ cm}^{-1}$ (valor guía), con sensor de temperatura Pt1000 integrado y con cable fijo para la conexión al 912/914 Meter.

Gracias al eje de plástico robusto e irrompible de PEEK, este sensor es altamente resistente a la tensión mecánica y es adecuado para medidas de conductividades intermedias (desde $15 \mu\text{S/cm}$ hasta 250 mS/cm), por ejemplo, en:

- Agua potable
- Aguas superficiales
- Aguas residuales

Accesorios opcionales

Order no.	Descripción
-----------	-------------

6.2301.060 KCl ion standard 250 mL

KCl ion standard ($c(\text{ion}) = 0.1000 \pm 0.0005 \text{ mol/L}$), can also be used as standard solution for conductometry 12.87 mS/cm (25°C).



6.2324.010 Conductivity standard $100 \mu\text{S/cm}$ 250 mL

Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$.



6.2324.110 Conductivity standard $100 \mu\text{S/cm}$, 5 x 30 mL

Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$.

