



Hamilton PRP-X300 - 250/4,0

6.1005.030

La columna de expulsión de iones Hamilton PRP-X300 es una columna de intercambio de cationes con una capacidad baja. La combinación de un copolímero de divinilbenceno-poliestireno con grupos de ácido sulfónico como intercambiador de iones es la solución ideal para los problemas de separación sencillos. Esta columna destaca por posibilitar la determinación de las sales de ácidos orgánicos, en especial la determinación de formiato muy sensible.

Accesorios opcionales

Order no.	Descripción
6.1011.030	Metrosep RP 2 Guard/3.5

The Metrosep RP 2 Guard/3.5 is a guard column for universal use. It reliably protects the analytical separation column against contaminations. It reliably protects the analytical separation column against contamination, removing the smallest particles, traces of iron oxide and bacteria. The Metrosep RP 2 Guard/3.5 helps to reduce costs; its filter disk can be replaced very easily.



6.1011.040

Metrosep RP 3 Guard HC/4.0

The Metrosep RP 3 Guard HC/4.0 is a guard column for universal use. It reliably protects the analytical separation column against contamination, securely removing lipophilic organic contamination, the smallest particles, traces of iron oxide and bacteria. The guard column is based on a polymer material and thanks to its larger pack volume, has a significantly higher capacity than the Metrosep RP 2Guard/3.5. It can be used throughout the entire pH range.



6.1011.130

Spare filter for RP 2 Guard/3.5, 10 pieces

