



Metrosep A Supp 5 - 150/4,0

6.1006.520

La versión de 150 mm de la Metrosep A Supp 5 destaca por sus excelentes propiedades de separación. Un alto número de platos y excelentes simetrías de picos facilitan los trabajos en la gama inferior de $\mu\text{g/L}$. El tamaño de partícula de 5 μm ayuda de manera decisiva a la eficacia de separación de esta columna. La Metrosep A Supp 5 - 150/4,0 ofrece la combinación óptima de selectividad y capacidad con la que se pueden solucionar de forma segura incluso complejas tareas de separación en poco tiempo. Esta propiedad hace que la Metrosep A Supp 5 - 150/4,0 sea una de las mejores columnas estándar de uso universal de la cromatografía iónica.

Accesorios opcionales

Order no.	Descripción
-----------	-------------

The Metrosep A Supp 5 Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contamination from the sample or the eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the respective separation column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 Guard/4.0 highly recommended.



The Metrosep A Supp 5 S-Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contaminations from the sample or eluent.

