



Metrosep A Supp 5 - 50/4,0

6.1006.550

En menos de 6 minutos, la Metrosep A Supp 5 - 50/4,0 separa los siete aniones estándar. Incluso el fluoruro está separado del pico de inyección y puede integrarse sin problemas. La columna, basada en un polímero de alcohol polivinílico, se caracteriza, como todas las columnas de la familia A Supp 5, por un alto número de platos y, en consecuencia, una excelente eficacia de separación. La Metrosep A Supp 5 - 50/4,0 es la columna adecuada cuando se deben solucionar tareas de separación sencillas en poco tiempo y sin renunciar a límites de detección muy bajos.

Accesorios opcionales

Order no.	Descripción
-----------	-------------

6.1006.500

Metrosep A Supp 5 Guard/4.0

The Metrosep A Supp 5 Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contamination from the sample or the eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the respective separation column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 Guard/4.0 highly recommended.



6.1006.540

Metrosep A Supp 5 S-Guard/4.0

The Metrosep A Supp 5 S-Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contaminations from the sample or eluent.

