



## Metrosep A Supp 7 - 150/2,0

6.1006.640

La Metrosep A Supp 7 - 150/2,0 es la más corta de las dos columnas A Supp 7 en la versión microbore. Con la Metrosep A Supp 7 - 150/2,0 es posible determinar estos iones hasta la gama de  $\mu\text{g/L}$  inferior de forma segura y precisa.

La alta sensibilidad de detección se logra usando el polímero de alcohol polivinílico de  $5\text{ }\mu\text{m}$  con el que pueden alcanzarse números de platos extremadamente elevados y con ellos excelentes propiedades de separación y detección. Además, es posible adaptar la separación a las exigencias específicas de la aplicación modificando la temperatura.

Esta columna microbore es especialmente apta para el uso con un detector MS.

### Accesorios opcionales

| Order no. | Descripción |
|-----------|-------------|
|-----------|-------------|

The Metrosep A Supp 5 Guard/2.0 reliably protects the 2 mm versions of the Metrosep A Supp 5 and Metrosep A Supp 7 IC anion columns against contamination from the sample or eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the respective separation column with virtually no dead volume ("On Column Guard System"). The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 Guard/2.0 highly recommended.



The Metrosep A Supp 5 S-Guard/2.0 reliably protects the 2 mm versions of the Metrosep A Supp 5 and Metrosep A Supp 7 IC anion columns against contamination from the sample or eluent.

It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is connected directly onto the respective separation column via capillary connection. The guard column prolongs the service life of the analytical column, with practically no influence on its chromatographic separating efficiency. The economical price and simple handling make using the A Supp 5 S-Guard/2.0 highly recommended.



The Metrosep A Supp 16 Guard/2.0 reliably protects the Metrosep A Supp 16 analytical separation columns with 2 mm inner diameter against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily onto the analytical column. No tools are required.



The Metrosep A Supp 16 S-Guard/2.0 reliably protects the Metrosep A Supp 16 analytical separation columns with 2 mm inner diameter against contamination. The guard column is connected to the separation column simply using capillary connections. No tools are required.

