



Metrosep C 3 - 100/4,0

6.1010.410

El novedoso material portador de alcohol polivinílico aumenta de forma significativa la capacidad de separación de cationes monovalentes y bivalentes. Las formas de los picos en esta columna de cationes son muy simétricas.

La columna de separación más corta de la familia Metrosep C 3 es particularmente apropiada para efectuar separaciones rápidas de los cationes estándar y separaciones de aminas orgánicas grandes.

Accesorios opcionales

Order no.	Descripción
6.1010.450	Metrosep C 3 Guard/4.0
	The Metrosep C 3 S-Guard/4.0 contains the C-3 column material and is used to protect Metrosep C-3 cation columns. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The Metrosep C 3 Guard/4.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.



The Metrosep C 3 S-Guard/4.0 contains the C-3 column material and is used to protect Metrosep C-3 cation columns. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column.

