



Metrosep A Supp 10 - 50/2,0

6.1020.250

La columna de separación Metrosep A Supp 10 - 50/2,0 está basada en un copolímero de poliestireno-divinilbenceno de alta capacidad, con un tamaño de partículas de solo 4,6 µm. Las propiedades de la columna se pueden adaptar específicamente a la aplicación actual mediante la temperatura, el flujo y la composición del eluyente. Las columnas de separación de 2 mm Metrosep A Supp 10 están hechas del mismo material que las columnas de separación de 4 mm correspondientes. La escasa longitud y la capacidad total relativamente baja de esta columna de 50 mm permite separaciones muy rápidas de aniones estándar.

La Metrosep A Supp 10 - 50/2,0 es ideal para sencillas tareas de separación y matrices no complicadas. Esta columna de separación microbore es adecuada para aplicaciones IC-MS gracias al reducido flujo.

Partes/accesorios 6.1020.250

Qt.	Order no.	Descripción
-----	-----------	-------------

2 PCS

6.2744.060

Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Accesorios opcionales

Order no.	Descripción
6.1020.600	<p>Metrosep A Supp 10 Guard/2.0</p> <p>The Metrosep A Supp 10 Guard/2.0 column reliably protects the A Supp 10 analytical microbore separation columns against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily and directly onto the analytical column without tools</p>

