



Metrosep A Supp 16 - 150/4,0

6.1031.420

La Metrosep A Supp 16 es ideal para problemas de separación de alta capacidad y posee una resolución excelente. La Metrosep A Supp 16 - 150/4,0 está basada en un copolímero de divinilbenceno-poliestireno con superficie funcionalizada. Los grupos funcionales están ligados de forma covalente.

La Metrosep A Supp 16 - 150/4,0 posee una resolución excelente y soluciona problemas de separación muy difíciles. Esta columna es ideal para aplicaciones que presentan una carga elevada de iones, pero que no requieren la máxima resolución. Es una de las columnas estándar en la cromatografía aniónica.

Partes/accesorios 6.1031.420

Qt.	Order no.	Descripción
-----	-----------	-------------

2 PCS

6.2744.060

Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Accesorios opcionales

Order no.	Descripción
6.1031.500	Metrosep A Supp 16 Guard/4.0

The Metrosep A Supp 16 Guard/4.0 reliably protects the Metrosep A Supp 16 analytical separation columns against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily onto the analytical column. No tools are required.



6.1031.510	Metrosep A Supp 16 S-Guard/4.0
------------	--------------------------------

The Metrosep A Supp 16 S-Guard/4.0 reliably protects the Metrosep A Supp 16 analytical separation columns against contamination. The guard column is connected to the separation column simply using capillary connections. No tools are required.

