



## Metrosep C 4 - 100/2,0

6.1050.210

La versión corta de la columna de separación de cationes microbore Metrosep C 4 de 2 mm de diámetro interior está concebida para determinaciones rápidas de los cationes estándar. Con ella se logran tiempos de retención muy cortos, a pesar de que los tiempos de elución del sodio y del amonio tienen una diferencia de 25 s. Si se usa un eluyente especial, con la Metrosep C 4 - 100/2,0 se pueden determinar los seis cationes litio, amonio, sodio, calcio, magnesio y potasio en menos de 5 minutos. Debido a su reducido flujo de eluyente, es especialmente apropiada para el acoplamiento IC-MS.

### Partes/accesorios 6.1050.210

| Qt. | Order no. | Descripción |
|-----|-----------|-------------|
|-----|-----------|-------------|

2 PCS

6.2744.060

Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



## Accesorios opcionales

| Order no.  | Descripción                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |                                                                                      |
|------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| 6.1050.600 | Metrosep C 4 Guard/2.0                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |   |
|            | <p>The Metrosep C 4 Guard/2.0 contains the C-4 column material and is used to protect all Metrosep cation columns with 2 mm inner diameters which have a substrate based on silica gel. Particles and contamination are reliably retained, considerably prolonging the service life of the analytical separation column. The economical price is an additional plus. The Metrosep C 4 Guard/2.0 also works based on the "On Column Guard System" and is screwed directly onto the separation column with virtually no dead volume.</p> |                                                                                      |
| 6.1050.610 | Metrosep C 4 S-Guard/2.0                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |  |