



Recipiente de reacción para medidas de estabilidad

6.1429.040

Recipiente de reacción para medidas de estabilidad con 743 Rancimat y 763 PVC Thermomat.
Embalaje con 117 unidades.

Accesorios opcionales


Order no.	Descripción
6.1808.090	Thread adapter M8 outer / M6 inner Outer thread M8, inner thread M6, e.g. for mounting 6.2418.XXX Air tubes in 6.1429.XXX Reaction vessels of stability measuring instruments.



6.2418.100

Air tube for oil measurements

Set of 117 pieces. In combination with 6.1429.040.



6.2418.120

Gas inlet tube for PVC measurements

Set of 117 pieces. In combination with 6.1429.040.

