



Cierre de rosca con septo, 100 uds.

6.1448.067

Cierre de rosca con septo, apto para frasco de muestras de 8 mL 6.2420.007, 100 uds.

Accesorios opcionales

Order no.	Descripción
6.1448.077	Septum, 100 pcs.

Replacement septum for 6.1448.067 screw cap, 100 pieces



6.2420.007

Sample glass with thread, 8 mL, 100 pieces

Sample glass with thread, 8 mL, suitable for 6.1448.067 screw cap, 100 pieces


