



Consola de soporte para 912/913/914

6.2001.130

Consola de soporte para covertir un 912/913/914 pH/Conductometer en un aparato de laboratorio.

Partes/accesorios 6.2001.130

Qt.	Order no.	Descripción
1 PCS	6.2013.010	Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2016.050

Support rod / 300 mm



1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors

