



Tapón septo

6.2730.020

Con junta tórica E.301.0041. Para septo 6.1448.010 y parte superior de recip. de titración 6.1414.030. Para aplicaciones KF.

Partes/accesorios 6.2730.020

Qt.	Order no.	Descripción
1 PCS	6.1448.010	Septum 12 mm, 5 pieces
		Set of 5 items.



1 PCS

6.2730.010

Screw nipple

Made of PP. For 6.1448.010 Septum or for mounting 6.1403.040 Drying tube or 6.0338.100 Electrode in 6.1414.030 and 6.1414.040 KF Titration vessel lids. Used for KF titrations.

