



Metrosep Amino Acids 1 - 100/4,0

6.4001.410

Columna de separación para la determinación de aminoácidos mediante detección UV/VIS tras derivatización post-columna con ninhidrina a 120 °C. La Metrosep Amino Acids 1 - 100/4,0 es la columna de separación estándar para los aminoácidos. La columna está basada en un material de poliestireno-divinilbenceno sulfonado. La determinación de los aminoácidos se efectúa mediante detección fotométrica tras una reacción post-columna con ninhidrina.

La Metrosep Amino Acids 1 - 100/4,0 permite separar hasta 44 aminoácidos, incluidos todos aquellos que se producen de forma natural en la investigación y en la práctica.

Partes/accesorios 6.4001.410

Qt.	Order no.	Descripción
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Accesorios opcionales

Order no.	Descripción	
6.1011.030	Metrosep RP 2 Guard/3.5	
	The Metrosep RP 2 Guard/3.5 is a guard column for universal use. It reliably protects the analytical separation column against contaminations. It reliably protects the analytical separation column against contamination, removing the smallest particles, traces of iron oxide and bacteria. The Metrosep RP 2 Guard/3.5 helps to reduce costs; its filter disk can be replaced very easily.	
6.1011.040	Metrosep RP 3 Guard HC/4.0	
	The Metrosep RP 3 Guard HC/4.0 is a guard column for universal use. It reliably protects the analytical separation column against contamination, securely removing lipophilic organic contamination, the smallest particles, traces of iron oxide and bacteria. The guard column is based on a polymer material and thanks to its larger pack volume, has a significantly higher capacity than the Metrosep RP 2Guard/3.5. It can be used throughout the entire pH range.	
6.1011.130	Spare filter for RP 2 Guard/3.5, 10 pieces	