



Juego de accesorios CI: diálisis inline Low Volume

6.5330.200

Set de accesorios para una diálisis inline rápida. Para el uso con el 858 Professional Sample Processor y una bomba peristáltica bicanal adicional.

Partes/accesorios 6.5330.200

Qt.	Order no.	Descripción
1 PCS	6.1602.150	Bottle attachment / GL 45 - 3 x UNF 10/32

For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.





Capillary for sample handling in IC.



For suppressor solutions and for acceptor solutions for inline dialysis and for inline ultrafiltration.



For the sample solution in inline dialysis.



1 PCS

6.1831.120

PEEK capillary 0.25 mm i.d., 45 cm



1 PCS

6.1831.160

PEEK capillary, 0.5 mm ID / 70 cm




1 PCS

6.2057.130

Dialysis cell holder

Holder for the dialysis cell



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



1 PCS

6.2714.010

Dialysis membrane (cellulose acetate)

Nominal pore size : 0.2 µm. Set of 50 pieces. For inline dialysis



1 PCS

6.2729.200

Low Volume dialysis cell

Including four 6.2744.060 stoppers. For rapid Inline Dialysis.



1 PCS

6.2744.000

Pressure screw PVDF

With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)



2 PCS

6.2744.034

Coupling olive/UNF 10/32, 2x

Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.



2 PCS

6.2744.070

Pressure screw short

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries




4 PCS

6.2744.160

Nozzle/UNF 10/32 coupling with security device

Security device for connecting pump tubing to a nozzle.



1 PCS

6.2755.000

Tubing cartridge

For the peristaltic pumps of IC instruments. Contact pressure adjustable.



1 PCS

6.2831.010

Plastic tweezers

For dialysis and filtration membranes

