



Juego de accesorios para la determinación de la corrección de temperatura en Biodiesel Rancimat.

6.5616.110

Set para el ajuste exacto de la temperatura

Partes/accesorios 6.5616.110

Qt.	Order no.	Descripción
1 PCS	6.1111.020	Pt100 temperature sensor for stability measuring devices

Pt100 temperature sensor (length 30.0 cm) with fixed cable and mini-DIN-plug. This glass-free temperature sensor is well suited for use with instruments for stability measurement (873 Biodiesel Rancimat and 893 Professional Biodiesel Rancimat).



2 PCS

6.2042.040

Clamp for temperature sensor

Holding clamp for Pt100 temperature sensor for stability measuring instruments



1 PCS

6.2326.000

Silicone oil for stability measuring instruments (50 mL)

Silicone oil for the determination of temperature correction in stability measuring instruments

