



OMNIS Coulometer sans agitateur

2.1018.0010

OMNIS Coulometer sans agitateur magnétique pour mode autonome ou comme pièce centrale d'un système de titrage OMNIS. Le coulomètre peut être complété par d'autres modules de titrage, de dosage ou de coulométrie.

La coulométrie est la méthode parfaitement adaptée à la détermination de la teneur en eau dans le domaine des traces (10 µg à 10 mg de teneur en eau absolue) pour les liquides, les matières solides et les gaz. En outre, le titrage coulométrique Karl Fischer est une méthode absolue, c'est pourquoi une détermination du titre est inutile

Contenu de la livraison 2.1018.0010

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

1.1018.0010

OMNIS Coulometer without stirrer

OMNIS Coulometer without magnetic stirrer



1 PCS

6.06001.970

Function license Coulometer

Function license for OMNIS Coulometer, includes

- Liquid Handling (LQH)
- Coulometric Karl Fischer titration (KFC)
- Coulometric determination of the bromine index according to ASTM D1492 (BRC)
- MEAS U/T/pH/CONC

MNI
MORE NEW LEVEL OF PERFORMANCE

Accessoires optionnels

Order no.	Description	
2.1009.0010	OMNIS Solvent Module OMNIS Solvent Module for reliable addition of solvent and aspirating off titration vessel content without contact with the reagents. New, improved version resistant against chloroform, chlorobenzene, xylene, and glacial acetic acid (requires OMNIS 4.4.1 or higher). Manual control or actuation via OMNIS Software for maximum user-friendliness. Including accessories, solvent, and waste bottle.	
6.02102.010	Cable MDL PL/SO 0.5 m Connecting cable for OMNIS instruments with MDL connector (0.5 m)	
6.02102.020	Cable MDL PL/SO 1 m Connecting cable for OMNIS instruments with MDL connector (1 m)	
6.02102.030	Cable MDL PL/SO 2 m Connecting cable for OMNIS instruments with MDL connector (2 m)	

6.02102.040 Cable MDL PL/SO 5 m

Connecting cable for OMNIS instruments with MDL connector (5 m)



6.02102.050 Cable MDL PL/SO 10 m

Connecting cable for OMNIS instruments with MDL connector (10 m)



6.05406.000 KFC equipment for OMNIS Coulometer with generator electrode with diaphragm

Complete accessory set for coulometric Karl Fischer titration.



6.05406.010 KFC equipment for OMNIS Coulometer with generator electrode without diaphragm

Complete accessory set for coulometric Karl Fischer titration.



6.06001.100 Function license DET

"Dynamic equivalence point titration" function license for one OMNIS titration system.



6.06001.110 Function license MET

"Monotonic equivalence point titration" function license for one OMNIS titration system.



6.06001.160 Function license KFT with conditioning

Function license for volumetric Karl Fischer titration with conditioning for an OMNIS titration system.



6.06001.320 Function license for 5 parallel titrations

Function license for a maximum of five parallel titrations using one OMNIS Titrator with additional Titration Modules or Dosing Modules.



6.1446.220

Reagent Replacement Solvent Module stopper

Stopper for reagent replacement in a coulometric titration cell with a Solvent Module.

