



OMNIS Coulometer avec agitateur

2.1018.0020

OMNIS Coulometer avec agitateur magnétique pour mode autonome ou comme pièce centrale d'un système de titrage OMNIS. Le coulomètre peut être complété par d'autres modules de titrage, de dosage ou de coulométrie.

La coulométrie est la méthode parfaitement adaptée à la détermination de la teneur en eau dans le domaine des traces (10 µg à 10 mg de teneur en eau absolue) pour les liquides, les matières solides et les gaz. En outre, le titrage coulométrique Karl Fischer est une méthode absolue, c'est pourquoi une détermination du titre est inutile.

Contenu de la livraison 2.1018.0020

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.06001.970

Function license Coulometer

Function license for OMNIS Coulometer, includes

- Liquid Handling (LQH)
- Coulometric Karl Fischer titration (KFC)
- Coulometric determination of the bromine index according to ASTM D1492 (BRC)
- MEAS U/T/pH/CONC



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2016.050

Support rod / 300 mm



Accessoires optionnels

Order no.	Description	
2.1009.0010	OMNIS Solvent Module	
	OMNIS Solvent Module for reliable addition of solvent and aspirating off titration vessel content without contact with the reagents. New, improved version resistant against chloroform, chlorobenzene, xylene, and glacial acetic acid (requires OMNIS 4.4.1 or higher). Manual control or actuation via OMNIS Software for maximum user-friendliness. Including accessories, solvent, and waste bottle.	
6.02102.010	Cable MDL PL/SO 0.5 m	
	Connecting cable for OMNIS instruments with MDL connector (0.5 m)	
6.02102.020	Cable MDL PL/SO 1 m	
	Connecting cable for OMNIS instruments with MDL connector (1 m)	
6.02102.030	Cable MDL PL/SO 2 m	
	Connecting cable for OMNIS instruments with MDL connector (2 m)	

6.02102.040 Cable MDL PL/SO 5 m

Connecting cable for OMNIS instruments with MDL connector (5 m)



6.02102.050 Cable MDL PL/SO 10 m

Connecting cable for OMNIS instruments with MDL connector (10 m)



6.05406.000 KFC equipment for OMNIS Coulometer with generator electrode with diaphragm

Complete accessory set for coulometric Karl Fischer titration.



6.05406.010 KFC equipment for OMNIS Coulometer with generator electrode without diaphragm

Complete accessory set for coulometric Karl Fischer titration.



6.06001.100 Function license DET

"Dynamic equivalence point titration" function license for one OMNIS titration system.



6.06001.110 Function license MET

"Monotonic equivalence point titration" function license for one OMNIS titration system.



6.06001.160 Function license KFT with conditioning

Function license for volumetric Karl Fischer titration with conditioning for an OMNIS titration system.



6.06001.320 Function license for 5 parallel titrations

Function license for a maximum of five parallel titrations using one OMNIS Titrator with additional Titration Modules or Dosing Modules.



6.1446.220

Reagent Replacement Solvent Module stopper

Stopper for reagent replacement in a coulometric titration cell with a Solvent Module.

