



786 Swing Head avec tête de transfert droite

2.786.0020

Accessoires supplémentaires pour le Robotic Sample Processor XL pour transférer / pipetter à partir de petits récipients échantillons dans des grands vases de titrage sur le portoir ou dans une cellule de titrage externe, pouvant être installée à droite de la station de travail.

Contenu de la livraison 2.786.0020

Qt.	Order no.	Description
1 PCS	1.786.0010	786 Swing Head with transfer head left Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.



1 PCS

6.1462.040

Robotic arm with transfer head for 786 Swing Head,
right swinging

Transfer robotic arm for liquid handling purposes



1 PCS

6.1562.100

Pipetting tubing / 10 mL

Pipetting tubing for fully automatic transfer of sample solutions.



1 PCS

6.1823.010

Guide shaft to pipetting tubing

For assembly in the transfer robotic arm on the Robotic Sample Processors



1 PCS

6.2621.120

Hexagon key 1.5 mm

