



786 Swing Head avec bras pivotant de titrage, vers la gauche ou vers la droite

2.786.0030

Accessoires supplémentaires pour le Robotic Sample Processor XL pour le titrage direct dans des récipients échantillons de 75 mL et plus.

Contenu de la livraison 2.786.0030

Qt.	Order no.	Description
1 PCS	1.786.0010	786 Swing Head with transfer head left
		Additional equipment for Robotic Sample Processors XL for transferring / pipetting from small sample vessels into larger titration vessels on the rack or into an external titration cell; this can be located to the left of the work station.



1 PCS

6.1462.050

Robotic arm with titration head for 786 Swing Head,
left or right swinging

Titration head for Swing Head



1 PCS

6.1808.040

Thread adapter M6 outer / M8 inner

Outer thread M6, inner thread M8.



1 PCS

6.1909.030

Stirring propeller / 95 mm

Propeller stirrer with PP shaft, PVDF tip. Length from the SGJ: 104 mm. For
KF applications with 722, 802 Rod Stirrer



2 PCS

6.2042.030

Buret tip clip (rod stirrer)

Clip for proper alignment of buret tips, dosing or aspiration tubings using
a rod stirrer in the titration beaker of a sample changer



