



823 Membrane Pump Unit – rinse

2.823.0030

Pompe à membrane pour la connexion directe aux Sample Processors. La version «Rincer» comprend de nombreux accessoires de rinçage.

Contenu de la livraison 2.823.0030

Qt.	Order no.	Description
1 PCS	1.823.0010	823 Membrane Pump Unit
		Membrane pump for direct connection to Sample Processors.



1 PCS

6.1621.000

Container 10 L

As rinsing or waste container in automated systems.



2 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



3 PCS

6.1805.060

FEP tubing / M6 / 60 cm

With light and kink protection



1 PCS

6.1808.170

Distributor for Sample Processors for rinsing and aspirating

Distributor for installation on Sample Processors for rinsing and aspirating when using externally connected pumps.



1 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



1 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



3 PCS

6.2740.020

Spray nozzle

For the fine-spraying of the rinsing solution.

