



Titrode iAg avec revêtement d'Ag₂S

6.00470.300S

Électrode annulaire en argent combinée, intelligente à membrane de verre pH comme électrode de référence et puce de mémoire intégrée pour les données du capteur.

L'anneau en argent est revêtu de sulfure (Ag₂S) pour une sensibilité accrue et une meilleure limite de sensibilité.



Cette électrode ne nécessitant aucune maintenance convient aux titrages par précipitation sous pH constant (nitrate d'argent comme titrant), par ex. de :

- Chlorure, bromure, iodure
- Sulfures
- Sulfure d'hydrogène
- Mercaptans
- Cyanures


Cette électrode se conserve dans de l'eau distillée.

Les iTrodes peuvent être utilisées sur un Titrand, un Ti-Touch ou un 913/914 Meter.

Contenu de la livraison 6.00470.300S

Qt.	Order no.	Description
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm Sleeve with SGJ 14/12 mm, polyethylene.
		
1 PCS	6.2008.040	Storage vessel Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.
		

Accessoires optionnels

Order no.	Description
2.854.0010	<div>854 iConnect</div> <div>The 854 iConnect – Electrode cable and measuring amplifier for intelligent "iTrodes" electrodes (cable length 150 cm).</div> <div></div>