



Cellule de mesure de conductivité $c = 0,5 \text{ cm}^{-1}$ avec Pt1000 (câble fixe, 0,65 m)

6.00923.080

Cellule de mesure de conductivité à 4 pôles avec constante de cellule $c = 0,5 \text{ cm}^{-1}$ (valeur guide) avec capteur de température intégré Pt1000 et câble fixe (0,65 m) de connexion à l'OMNIS Measuring Module Conductivity.

Grâce à la tige en PEEK robuste/incassable, ce capteur présente une très grande résistance mécanique et est adapté aux mesures de moyenne conductivité (de $15 \mu\text{S/cm}$ à 250 mS/cm) comme dans :

- l'eau potable
- les eaux de surface
- les eaux usées

Contenu de la livraison 6.00923.080

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



Accessoires optionnels

Order no.	Description	
6.02101.020	Measuring module conductivity	
	Measurement channel for one OMNIS Titrator or Titration Module for the connection of conductivity measuring cells.	
6.2301.060	KCl ion standard 250 mL	
	KCl ion standard ($c(\text{ion}) = 0.1000 \pm 0.0005 \text{ mol/L}$), can also be used as standard solution for conductometry 12.87 mS/cm (25°C).	
6.2324.010	Conductivity standard $100 \mu\text{S/cm}$ 250 mL	
	Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$.	
6.2324.110	Conductivity standard $100 \mu\text{S/cm}$, 5 x 30 mL	
	Conductivity standard for calibration of conductivity measuring cells with cell constant = $0.1/\text{cm}$.	