



Metrosep A Supp 17 - 100/4,0

6.01032.410

Les colonnes Metrosep A Supp 17 sont des colonnes de séparation des anions pour un emploi à température ambiante, présentant un rapport performance/prix très avantageux. La Metrosep A Supp 17 - 100/4,0 permet de séparer les anions standard dans des délais brefs. Grâce à la grande flexibilité des débits d'écoulement (jusqu'à 1,8 mL/min), il est possible d'obtenir des temps d'analyse très réduits, en fonction des problèmes de séparation.

Contenu de la livraison 6.01032.410

Qt.	Order no.	Description
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Accessoires optionnels

Order no.	Description	
6.01032.500	Metrosep A Supp 17 Guard/4.0 The Metrosep A Supp 17 Guard/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily onto the analytical column. No tools are required.	
6.01032.510	Metrosep A Supp 17 S-Guard/4.0 The Metrosep A Supp 17 S-Guard/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. The guard column is simply connected to the separation column using capillary connections. No tools are required.	
6.01032.530	Metrosep A Supp 17 S-Guard - 50/4.0 The Metrosep A Supp 17 S-Guard 50/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. The guard column is simply connected to the separation column using capillary connections. No tools are required.	