



Metrosep C 6 - 250/2,0

6.01051.230

La colonne Metrosep C 6 - 250/2,0 est la colonne microbore pour cations avec la plus haute capacité de la série C 6 . Elle est destinée aux applications exigeant des performances de séparation très élevées. Cette colonne permet d'analyser de manière fiable des échantillons présentant des différences de concentration extrêmes. La séparation en particulier du sodium et de l'ammonium est excellente.

La colonne convient pour le couplage CI-SM.

Contenu de la livraison 6.01051.230

Qt.	Order no.	Description
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Accessoires optionnels

Order no.	Description
6.01051.600	Metrosep C 6 Guard/2.0

The Metrosep C 6 Guard/2.0 contains the C 6 column material and is used to protect against particles and contamination. This considerably lengthens the service life of the analytical separation column. The Metrosep C 6 Guard/2.0 works according to the "On Column Guard System" and is screwed directly onto the respective separation column with virtually no dead volume.

