



Unité de cylindre OMNIS 20 mL

6.03001.220

Unité de cylindre intelligente 20 mL pour un OMNIS Titrator, Titration Module ou Dosing Module. Tuyaux de dosage et pointe anti-diffusion compris.

Contenu de la livraison 6.03001.220

Qt.	Order no.	Description
1 PCS	6.01503.220	OMNIS 20 mL cylinder unit without accessories

Intelligent 20 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.



1 PCS

6.1446.030

Ball stopper / B-14/(15)

For the movable holding of buret tips in titration vessels



2 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1543.200

Titration tip / M6 thread

Titration tip with antidiffusion valve.



2 PCS

6.1805.100

FEP tubing / M6 / 40 cm

With light and kink protection



1 PCS

6.2008.030

Buret tip storage vessel for 807 Dosing Unit

Together with the 6.2008.050 storage vessel holder, it provides a holder for the buret tip on the 807 Dosing Unit.



1 PCS

6.2709.070

Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders



1 PCS

6.2739.000

Wrench

For tightening connectors



Accessoires optionnels

Order no.	Description	
6.05700.140	Consumable Kit, cylinder unit OMNIS	
	FEP tubing connections M6 (40 cm) and antidiffusion tip for 6.03001.xxx and 6.03004.xxx OMNIS cylinder units.	
6.05700.431	Consumable Kit, cylinder unit OMNIS 20 mL	
	Cylinder unit OMNIS 20 mL complete with cylinder, piston, cylinder base, valve disk, centering tube, and cylinder housing, including programmed data chip.	