



Équipement CI : boucle complète Full-Loop, PP, 889 Basic

6.05330.220

L'équipement CI comprend tous les accessoires nécessaires au transfert d'échantillons du 889 IC Sample Center Basic par pompe péristaltique au chromatographe ionique.

Contenu de la livraison 6.05330.220

Qt.	Order no.	Description
1 PCS	6.1602.150	Bottle attachment / GL 45 - 3 x UNF 10/32

For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.



1 PCS

6.1608.070

Eluent bottle / 2 L / GL 45



1 PCS

6.1803.030

PTFE capillary 0.5 mm i.d. / 3 m

Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor,
IC Liquid Handling Dialysis Unit



1 PCS

6.1826.360

Pump tubing LFL (white/white), 3 stoppers

For Sample Processors



1 PCS

6.1841.000

PEEK transfer capillary 2 mL, 5 m



1 PCS

6.2744.034

Coupling olive/UNF 10/32, 2x

Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.



1 PCS

6.2744.070

Pressure screw short

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries

